


 **Bottle Top Dispenser FORTUNA UNIVERSAL® OPTIFIX® HF** made of PTFE, with PTFE-encased dispensing piston and PTFE-cylinder, FIX adaptor made of solid PTFE, comes with 2 PTFE adaptors GL 40+45, 1 PTFE intake tube, single pack

 **Flaschenaufsatzdispenser FORTUNA UNIVERSAL® OPTIFIX® HF** aus PTFE, mit PTFE-ummanteltem Dosierkolben und PTFE-Dosierzylinder, FIX-Adapter aus massivem PTFE, inkl. 2 PTFE-Adapter GL 40+45, 1 PTFE-Ansaugschlauch, VE 1

 **NEW**  
Not spoilt for choice any longer! One dispenser for (nearly) all liquids!

 **NEU**  
Keine Qual der Wahl mehr! Einer für (fast) alle Flüssigkeiten!

| Volume<br>Volumen | grad.<br>grad. | Art.-No.<br>Art.-Nr. |
|-------------------|----------------|----------------------|
| 1 - 5 ml          | 0.1 ml         | 101.106-33           |
| 2 - 10 ml         | 0.2 ml         | 101.106-38           |
| 6 - 30 ml         | 0.5 ml         | 101.106-44           |
| 10 - 50 ml        | 1.0 ml         | 101.106-48           |

suitable bottles on page 233 · passende Flaschen auf Seite 233

 For all liquids that can be handled with the OPTIFIX dispensers and for hydrofluoric acid

 Für alle mit den OPTIFIX Dispensern dosierbaren Medien und für Flusssäure