


 **Bottle Top Dispenser**
FORTUNA® OPTIFIX® SAFETY,
 with PTFE-encased dispensing
 piston with sealing lip,
 FIX adaptor made of solid PTFE,
 comes with 2 PTFE adaptors
 GL 40+45, 1 PTFE intake tube,
 single pack

 **Flaschenaufsatzdispenser**
FORTUNA® OPTIFIX®
SAFETY, mit
 PTFE-ummanteltem
 Dosierkolben mit Dichtlippe,
 FIX-Adapter aus massivem
 PTFE, inkl. 2 PTFE-Adapter
 GL 40+45, 1 PTFE-Ansaug-
 schlauch, VE 1



ISO 8655

Volume	grad.	Accuracy	Tolerance	Art.-No.
Volumen	grad.	Richtigkeit	Toleranz	Art.-Nr.
0.4 - 2 ml	0.1 ml	$\leq \pm 0.6 \%$	$\leq \pm 0.2 \%$	101.107-27
1 - 5 ml	0.1 ml	$\leq \pm 0.6 \%$	$\leq \pm 0.2 \%$	101.107-33
2 - 10 ml	0.2 ml	$\leq \pm 0.6 \%$	$\leq \pm 0.2 \%$	101.107-38
6 - 30 ml	0.5 ml	$\leq \pm 0.6 \%$	$\leq \pm 0.2 \%$	101.107-44
10 - 50 ml	1.0 ml	$\leq \pm 0.6 \%$	$\leq \pm 0.2 \%$	101.107-48
20 - 100 ml	2.0 ml	$\leq \pm 0.6 \%$	$\leq \pm 0.2 \%$	101.107-51

 For dispensing strong acids and aggressive chemicals
 Zur Dosierung von starken Säuren und aggressiven Chemikalien